Air Mass Flow Sensor

Release 0.0.1

GianAndrea Mueller

Apr 04, 2021

CONTENTS

1	Troubleshooting	3
2	Acknowledgements 2.1 Installation 2.2 Setup	5 . 5 . 6
	2.3 Drivers	. 10 . 15 . 17
Ру	thon Module Index	19
In	lex	21

This project contains the software for a student experiment at ETH Zürich based on sponsorship by Sensirion AG. The student experiment aims to illustrate the sensor principle of a thermal air mass flow meter, which is investigated mainly from the control point of view. The students are required to implement control algorithms for the regulation of the heater.

The software provided here mainly does two things:

- 1. Communication with all attached devices:
 - 1. Sensirion SHT temperature sensors connected to a Sensirion Sensor Bridge
 - 2. A Sensirion SFM massflow meter
 - 3. A custom built heater being driven with a PWM signal

This task is taken over by the *Setup* class. It handles all interactions with the hardware and for this purpose makes use of the different drivers, as seen on page *Drivers*.

- 2. Allowing interactions:
 - 1. Displaying the current system status
 - 2. Walking the student through different steps of the experimentation
 - 3. Handling interactions with the setup

These tasks are solved with a PyQt5 based graphical user interface as described in section GUI.

CHAPTER

ONE

TROUBLESHOOTING

When experiencing issues with soft- or hardware consider section FAQ.

CHAPTER

TWO

ACKNOWLEDGEMENTS

The icons in the GUI are made by Yusuke Kamiyamane and used under CC BY 3.0.

The GUI frontend, QT 5.0 is used under LGPL 3.0.

2.1 Installation

2.1.1 Windows 10

Setup GUI

- 1. Install python 3.8 or newer.
- 2. Set your PATH variable such that it includes the Scripts folder of your python installation.
- 3. Go to 01_SETUP/WINDOWS and run py -m setup in the cmd shell.
- 4. Install Sensirion Control Center to allow the sensor bridge to communication with the computer. Important: Select yes when asked for driver installation at the end of the process.
- 5. Find the finished executable at 02_SOFTWARE/disp.

Setup Sensirion USB Sensor Viewer

- 1. Install the Sensirion USB Sensor Viewer.
- 2. Select COM HARDWARE: RS485/USB Sensor Cable.
- 3. Select Sensor Product: DP Sensors (SDP3x/SDP8xx).
- 4. Execute *Drivers/identify_differential_pressure_sensor.py* with a local python environment. This will give you an overview of all connected sensors and print the comport ID of the pressure sensor in the final line.
- 5. Enter the previously found comport number in the RS485 Sensor viewer and connect.

Debugging

The installation is based on pyinstaller. It is configured via the 02_SOFTWARE/main.spec file. Set debug=True and console=True to receive informative output on the cmd shell upon launching the program.

2.2 Setup

class setup.**Setup** (*config: Utility.ConfigurationHandler.ConfigurationHandler*) The Setup handles all interaction with the hardware of the experiment.

Parameters

- serials (dict) Dictionary of device names and corresponding USB serials.
- **t_sampling_s** (*float*) Measurement sampling time in seconds.
- **interval_s** (*float*) Total buffered time interval in seconds, which in combination with the sampling time defines the number of stored measurements.

 $_measure_normal_mode() \rightarrow dict$

Measures all devices.

Returns A dictionary with all measured signals.

_measure_simulation_mode() \rightarrow dict

When no devices are connected random values are generated instead of actual measurements.

Returns A dictionary with all signals

 $_setup_measurement_buffer() \rightarrow Utility.MeasurementBuffer.MeasurementBuffer$

Defines the set of recorded signals and creates a corresponding MeasurementBuffer.

Returns An instance of MeasurementBuffer containing a deque instance for every signal.

See also:

Module Utility.MeasurementBuffer.MeasurementBuffer

$\texttt{close()} \rightarrow None$

Closes all connected devices.

 $\texttt{disable_output()} \rightarrow None$

Disable the output for either pwm or pid mode.

enable_output (*desired_pwm_output=0*) \rightarrow None

Enables the output in either pwn or pid mode.

Parameters desired_pwm_output (*float*) – Optionally enable pwm mode with a predefined nonzero output.

get_current_flow_value()

Getter for last set target flow value.

Returns Last set target flow value in normalized units.

$\texttt{measure()} \rightarrow None$

Handles measuring and storing signals depending on the system mode and handles updating the PID controller output.

See also:

```
_measure_simulation_mode() _measure_normal_mode() Utility.
MeasurementBuffer.MeasurementBuffer
```

open () \rightarrow None

Finds and opens all the USB devices previously defined within *self.serials* by their serial number. If one of the devices is not responsive or cannot be found, the setup is switching to simulation mode in which all measurements are simulated. This allows to test the GUI without any attached devices.

See also:

Module Drivers.DeviceIdentifier.DeviceIdentifier

```
\texttt{reset\_temperature\_calibration()} \rightarrow None
```

Reset the current temperature offset to zero.

```
reverse\_temp\_sensors(update=True) \rightarrow None
```

Reverse the order of the temperature sensors if the have been set up wrongly.

save_measurement_buffer (folder, name, type='mat')

Saves the current measurement buffer to a file. :param folder: Destination folder. :param name: Name of the file. A time tag will be appended for uniqueness. :param type: To allow different export filetypes.

set_flow(value)

Interface to the SFC5xxx drive for defining the current flow setpoint

Parameters flow (*float*) – The desired massflow in normalized units, in [0, 1].

$set_kd(kd: float) \rightarrow None$

Allows setting the Kd gain of the controller.

Parameters kd (float) - Kd gain of the controller

set_ki (*ki: float*) \rightarrow None

Allows setting the Ki gain of the controller.

Parameters ki (float) - Ki gain of the controller

set_kp (kp: float) \rightarrow None

Allows setting the Kp gain of the controller.

Parameters kp(float) - Kp gain of the controller.

set_pid_parameters (kp=None, ki=None, kd=None) \rightarrow None Interface to the pid-setting functionality of simple_pid.

Parameters

- **kp** (float) Kp gain of the controller
- ki (float) Ki gain of the controller
- kd (float) Kd gain of the controller

```
set_pwm(value: float) \rightarrow None
```

Safely sets the desired PWM value depending on the current system mode.

Parameters value (*float*) – Desired PWM value as a normalized value between 0 and 1.

See also:

setup.Mode

set_setpoint (*value: float*) \rightarrow None

Allows to define the temperature difference setpoint.

Parameters value (*float*) – Positive value smaller 20 degrees.

$\texttt{set_temperature_calibration()} \rightarrow None$

Record the current temperature offset, assuming steady state.

 $start_buffering() \rightarrow None$

Start recording measurements in the MeasurementBuffer and delete previously recorded measurements.

```
start_direct_power_setting () \rightarrow None Start pwm mode with the output set to off.
```

```
start measurement thread() \rightarrow None
```

Creates a thread.Timer that schedules future measurements at the desired sampling time.

See also:

Utility.Timer.RepeatTimer

start_pid_controller (*setpoint=None*) \rightarrow None Start pid mode with the output set to off.

Parameters setpoint (*float*) – Can be used to define a new temperature difference setpoint.

stop_buffering () \rightarrow None Stop recording measurements in the MeasurementBuffer.

2.2.1 Setup Modes

class setup.Mode(value)

Defines a set of system modes.

- 1. IDLE: Before any of the experiment modes has been loaded the system is idle.
- 2. FORCE_PWM_OFF: In this mode the pwm can be set directly, but the output is currently turned off.
- 3. FORCE_PWM_ON: In this mode the pwm can be set directly.
- 4. PID_OFF: In this mode the pid parameters can be set, but the output is currently turned off.
- 5. PID_ON: In this mode the pid parameters can be set and the controller is allowed to set pwm values.

2.2.2 Additional Utility

Logging Facility

Utility.Logger.setup_custom_logger(name: str, level: int) \rightarrow logging.log Sets the logging format, level and name of the logger.

Parameters

- **name** (*str*) Name of the logger.
- **level** (*int*) Initial logging level.

Returns Returns a log.

Measurement Buffer

The MeasurementBuffer holds a number of deque instances, one for each recorded signal and manages them as a ring buffer, always keeping a record of the most up to date measurements, reaching back *buffer_interval_s* seconds.

Parameters

- **signals** (*list*) List of signal names.
- **sampling_time_s** (*float*) Measurement sampling time in seconds.
- **buffer_interval_s** (*float*) Total buffered time interval in seconds which together with the sampling time defines the number of measurements to be stored.

$\texttt{clear()} \rightarrow None$

Clears the buffer.

update (*measurement: dict*) \rightarrow None

A buffer update is done by adding an entry to each signal buffer. Before the buffer is full this leads to an increase in length, afterwards the deque instances automatically forget their oldest entry in favor of the new one.

Parameters measurement (*object*) – Dictionary containing a value for each signal name

Timer

class Utility.Timer.RepeatTimer(interval, function, args=None, kwargs=None)

The RepeatTimer is a special type of timer thread that can be run indefinitely and executes a given function each time a specified interval has passed.

Note: Example of usage:

```
def dummyfn(msg="foo"):
    print(msg)
timer = RepeatTimer(interval=1, function=dummyfn)
timer.start()
time.sleep(5)  # During which 5 calls of dummyfn will happen.
timer.cancel()
```

$\texttt{run} \text{ () } \to \text{None}$

Method representing the thread's activity.

Overrides *Timer.run* such that we have a repeated timer.

2.3 Drivers

2.3.1 Device Identifier

class Drivers.DeviceIdentifier.DeviceIdentifier(serials: dict)

The DeviceIdentifier lists all connected USB devices and tries to identify all devices listed in *self.serials* with their respective serial port, which are subsequently available as *self.serial_ports*

Parameters serials (dict) – Dictionary with USB names as keys and USB serials as values.

Note: If the USB serials are unknown when launching the program first simply supply a dictionary with placeholders. DeviceIdentifier supplies information on all available devices upon failing to find one of the devices in the serials dictionary.

Warning: Windows detects USB serials differently than Linux. As experienced in the creating of this software, a serial read on a Linux system must be appended with the letter 'A' to be detected on a Windows system. To offer platform independence the serials must be given in 'Linux-Form' and are automatically appended with the letter 'A' when the program is run on Windows.

open()

Detects the current os. For Windows the letter 'A' is appended to the Linux-specific serial of the device. For Linux an additional tty-setup script is executed to allow detection of all USB devices. After that the serials of the available devices are compared and linked to *self.serials*.

Returns Returns True if all devices listed in self.serials could be found, False otherwise.

2.3.2 Platforms

Platform Base

class Drivers.PlatformBase.**PlatformBase** (*name: str*) Abstract base class for all platforms used in this project.

Parameters name (*str*) – Each platform must have a unique name.

 $\texttt{close()} \rightarrow None$

Disconnects the platform if it is currently connected.

 $\textbf{open}~(~)~\rightarrow bool$

Attempts to connect the platform and reports success.

Returns True if connected succesifully, False otherwise.

Shdic IO Module - The Heater

The ShdlcIoModule represents the custom Sensirion HDLC IO Box that allows driving the heater with a PWM output.

Parameters

- **serial_port** (*str*) Comport the IO box is connected to
- **baudrate** (*int*) Baudrate of the connection
- slave_address (int) Slave address
- input_pins (list) list of integers of the input pins

```
connect () \rightarrow bool
```

Attempts to connect the ShdlcIoModule

Returns True if connected successifully, False otherwise.

disconnect () \rightarrow None Sets all outputs off.

get_analog_input () \rightarrow float Measures actual voltage on ADC input

Returns A voltage between 0-10V

 $\begin{array}{l} \texttt{get_analog_output}() \rightarrow \texttt{float}\\ \texttt{Gets actual voltage for DAC output} \end{array}$

Returns A voltage between 0-10V

get_digital_io (*io_bit: int*) \rightarrow bool Reads a digital io pin.

Parameters io_bit (*int*) – Output bit to read

Returns True if digital bit is set.

get_pwm (pwm_bit: int) \rightarrow int Reads the current pwm setting.

Parameters pwm_bit (int) – The index of the PWM channel to be used (0, 1)

Returns A duty cycle value between 0 - 65535

$\texttt{is_connected()} \rightarrow bool$

Attempts to read the serial number of the device to check if it is connected.

Returns True if connected, False otherwise.

- $$\label{eq:set_all_digital_io_off()} \begin{split} & \texttt{set_all_digital_io_off()} \to None \\ & \text{Turns of all digital pins.} \end{split}$$
- **set_analog_output** (*value: int*) \rightarrow None Sets the analog output

Parameters value (*int*) – A voltage between 0-10V

set_digital_io (*io_bit: int, value: bool*) \rightarrow None Sets a digital output pin.

Parameters

- io_bit (int) The digital pin index to set
- value (bool) True if set to on

set_pwm (*pwm_bit: int, dc: int*) \rightarrow None Sets the pwm output

Parameters

• pwm_bit (*int*) – The index of the PWM channel to be used (0, 1)

• dc (int) -

Returns A duty cycle value between 0 - 65535

Sensirion Sensor Bridge (EKS)

```
class Drivers.SHT.EKS(serial_port: str)
```

Bases: Drivers.PlatformBase.PlatformBase

EKS represents a Sensirion Sensor Bridge which is used to communicate to a range of sensor via I2C.

Parameters serial_port (*str*) – Name of the port to which the EKS is connected.

```
\begin{array}{l} \textbf{connect} (\,) \ \rightarrow bool \\ Attempts to connect to the EKS. \end{array}
```

Returns True if connected sucessifully, otherwise the encountered exception will be returned.

```
connect_sensors ( ) \rightarrow None Attempts to connect sensors at both EKS ports.
```

- disconnect () \rightarrow None Closes all connected sensors.
- is_connected () \rightarrow bool Tests if the EKS is responsive.

Returns True if the EKs serial number can be read, False otherwise.

```
\texttt{measure()} \rightarrow list
```

Measures both channels if a sensor is attached

Returns A list of measured values.

2.3.3 Sensors

Sensor Base

```
class Drivers.SensorBase.SensorBase(name)
```

Abstract base class for all sensors used in this project.

Parameters name (*str*) – Each sensor must have a unique name.

 $\texttt{close()} \rightarrow None$

Disconnects the sensor if it is currently connected.

 $open() \rightarrow bool$

Attempts to connect the sensor and reports success. :return: True fi connected successifully, False otherwise

Sensirion Humidity Temperature (STH)

class Drivers.SHT.**SHT** (*device_port: int, shdlc_device: sensirion_shdlc_sensorbridge.device.SensorBridgeShdlcDevice, name='SHT'*)

Bases: Drivers.SensorBase.SensorBase

SHT represents either an SHT85 or an STH31 of the Sensirion Humidity Temperature (SHT) sensor range, connected via the Sensirion Sensor Bridge (EKS).

Parameters

- device_port (SensorBridgePort) EKS port, either ONE or TWO.
- **shdlc_device** (SensorBridgeShdlcDevice) Instance of the controlling EKS.
- **name** (*str*) Name of the sensor.

_convert_humidity (*data: bytearray*) \rightarrow float

Converts the raw sensor data to the actual measured humidity according to the data sheet Sensirion_Humidity_Sensors_SHT3x

Parameters data (*bytearray*) – 2 bytes, namely number 4 (humidity MSB) and 5 (humidity LSB) of the answer delivered by the sensor.

Returns The relative humidity measured by the sensor in percent.

__convert_temperature (*data: bytearray*) \rightarrow float

Converts the raw sensor data to the actual measured temperature according to the data sheet Sensirion_Humidity_Sensors_SHT3x

Parameters data (*bytearray*) – 2 bytes, namely number 1 (temperature MSB) and 2 (humidity LSB) of the answer delivered by the sensor.

Returns The temperature measured by the sensor in degrees Celsius.

$\texttt{connect} \text{ () } \rightarrow \text{bool}$

Attempts to connect the sensor and signals success by blinking the corresponding port's LEDs.

Returns Returns True if connected successifully, False otherwise.

connect_sensor (*supply_voltage: float, frequency: int*) \rightarrow None

Connection of a sensor attached to the sensirion sensor bridge according to the quick start guide to sensirion-shdlc-sensorbridge.

- Parameters
 - **supply_voltage** (*float*) Desired supply voltage in Volts.
 - **frequency** (*int*) I2C frequency in Hz

$\texttt{disconnect}() \rightarrow None$

Called by SensorBase.close upon deletion of this class. Switches supply off.

$\texttt{is_connected()} \rightarrow bool$

Check if the sensor operates correctly

Returns True if the status register can be read, False otherwise

$\texttt{measure()} \rightarrow dict$

Implementats a single shot measurement according to the SHT3x datasheet. A high repeatability measurement with clock stretching enabled is performed.

Returns Dictionary containing temperature in degrees Celsius and relative humiditiy in percent.

 $\begin{array}{l} \textbf{read_status_reg()} \rightarrow by tearray \\ Reads the status register \end{array}$

Returns Status register value as bytearray.

Sensirion Mass Flow Meter / Controller (SFM / SFC)

```
class Drivers.SFX5400.SFX5400(serial_port: str, name='Sfc5400')
Bases: Drivers.SensorBase.SensorBase
```

SFX5400 represents either a Sensirion Flow Controller (SFC) or a Sensirion Flow Meter (SFM) of type 5400.

Parameters

- **serial_port** (*str*) Name of the comport the SFX is connected to.
- **name** (*str*) Name of the device.

$\texttt{connect}() \rightarrow \texttt{bool}$

Attempts to connect to the SFX and reports success.

Returns True if connected successifully, False otherwise.

$\texttt{disconnect}() \rightarrow None$

Disconnects the device.

get_device_information (*index: int*) \rightarrow str

Retrieves device information depending on the index given.

Parameters index (*int*) – Integer between 1 and 3 to request on of the data below:

- 1. Product Name
- 2. Article Code
- 3. Serial number

Returns String containing the requested information.

is_connected()

Checks if the device is connected by reading its serial number.

Returns True if connected, False if not.

$\texttt{measure()} \rightarrow dict$

Measures the current mass flow.

Returns Dictionary containing the measurement.

$\texttt{set_flow} (\textit{setpoint_normalized: float}) \rightarrow \texttt{bool}$

Sets the current desired mass flow if a flow controller is connected.

Parameters setpoint_normalized (float) – Flow setpoint as normalized input between 0 and 1.

Returns True if set successifully, False if exception occured.

2.4 GUI

The structure of the graphical user interface can be described as follows: The outermost layer is within the main file, which deploys the Qt application and loads the main window.

The main window then controls the different experiment pages GUI.ExperimentPages.ExperimentPage, one for each experimentation step, with a stacked layout and manages the switching between those pages. The pages are built up from a series of widgets as defined in sections *Widgets* and *Live Plots*.

2.4.1 Main Window

```
class GUI.MainWindow.MainWindow(setup: setup.Setup, *args, **kwargs)
Defines the main window of the application.
```

Parameters setup (Setup) – Instance of Setup to allow access to sensors and actuators.

_calibrate_temperature() \rightarrow None

Toolbar action; Allows to set the current delta T to zero by saving the current temperature difference and subtracting it from the second temperature measurement.

_change_competition_mode() \rightarrow None

Toolbar action; Allows to set the current view to competition mode.

_go_to_next_view() ightarrow None

Toolbar action; Switches to the next view in the main layout stack.

_go_to_previous_view() ightarrow None

Toolbar action; Switches to the previous view in the main layout stack.

 $_$ **reset_plots**() \rightarrow None

Toolbar action; Allows to reset all visible plots to their original view.

_reset_temperature_calibration() \rightarrow None

Toolbar action; Allows to reset the calibration temperature difference to zero.

$_\texttt{reverse_temperature_sensors()} \rightarrow None$

Menu action; Allows to switch the order of the temperature sensors if the hardware setup is the wrong way around.

_save_measurement_buffer()

Toolbar aciton; Allows to save the measurement buffer as a Matlab .mat file.

 $_$ start $_$ recording() \rightarrow None

Toolbar action; Allows to restart recording measurements. Clears the buffer.

```
\_stop\_recording() \rightarrow None
```

Toolbar action; Allows to stop recording measurements and thus freeze the plots. :return:

 $_toggle_massflow(state=None) \rightarrow None$

Toolbar action; Allows to turn the massflow output on or off

Parameters state (bool) - Set True to turn the output state to on, or False vice versa.

_toggle_output (*state=None*) \rightarrow None

Toolbar action; Allows to turn the pwm output on or off.

Parameters state (bool) - Set True to turn the output state to on, or False vice versa.

_toggle_setpoint()

Toolbar action; Allows to change the temperature difference setpoint

 $setup_status_bar() \rightarrow None$

Sets up a status bar displaying sponsor layouts and tips for hovered over widgets.

```
\texttt{setup\_tool\_bar()} \rightarrow None
```

Adds a toolbar to the main window and defines a set of actions for it.

2.4.2 Widgets

class GUI.CustomWidgets.Widgets.**FancyPointCounter** (*setup*, **args*, ***kwargs*) Bases: PyQt5.QtWidgets.QLCDNumber

Custom version of the QLCDNumber.

property value

class	GUI.CustomWidgets.Widgets.CompetitionWidget	(setup:	setup.Setup,
		start_recording_action	: Callable,
		stop_recording_action.	: Callable,
		enable_output_action:	Callable,
		*args, **kwargs)	

Bases: GUI.CustomWidgets.BaseWidgets.FramedWidget

The CompetitionWidget allows to start a recording of the current performance and displays the number of points reached.

_update_process_values (*running_time_s*) \rightarrow None Container function for updates that are specific to the inheriting widgets

 $_update_progress() \rightarrow None$

Update the progressbar to show the current remaining time. If the recording interval has passed the process is stopped.

```
reset () \rightarrow None
Reset the competition widget upon reloading it.
```

```
class GUI.CustomWidgets.Widgets.StatusWidget(setup, *args, **kwargs)
    Bases: GUI.CustomWidgets.BaseWidgets.FramedWidget
```

The StatusWidget displays the current temperatures, flow and temperature difference measured.

```
_update_lcds () \rightarrow None
Updates the displayed values.
```

2.4.3 Live Plots

```
class GUI.CustomWidgets.LivePlots.LivePlotSignal (name: str, identifier: str, color: str,
width=1)
```

A LivePlotSignal stores all the information needed to identify and plot a single signal.

Parameters

- **name** (str) Name of the signal, to be displayed on the legend of the plot the signal is shown on
- **identifier** (str) Identifier of the signal, used to retrieve the signal from the measurement buffer of the setup
- color (str) Color of the plotted line used to instantiate the corresponding pen
- width (float) Width of the plotted line used to instantiate the corresponding pen

Note: Selecting integer values for the width parameter results in smoother plots.

The LivePlotWidget makes use of pyqtgraph to allow plotting a number of signals. It automatically updates.

Parameters

- setup (Setup) Instance of the current setup to allow access to the measurement buffer
- title (*str*) Title of the plot
- ylabel (str) Label of the y-axis
- ylims (Tuple) Limits of the y-axis

add_signals (signals: list) \rightarrow None Add a list of signals to the plot.

Parameters signals (*list*) – List of LivePlotSignals

update_plot_data() Handles the updating of a LivePLotWidget.

```
class GUI.CustomWidgets.LivePlots.LivePlotWidgetCompetition(setup: setup.Setup,
```

title, *ylabel*, *ylims*, **args*, ***kwargs*)

Bases: GUI. CustomWidgets. LivePlots. LivePlotWidget

Specialized LivePLotWidget allowing only two signals and adding color between the two corresponding lines. Used to visualize the integral of the control error.

add_signals (*reference_signal*, *actual_signal*) Add a list of signals to the plot.

Parameters signals (*list*) – List of LivePlotSignals

update_plot_data()

Handles the updating of a LivePLotWidget.

```
class GUI.CustomWidgets.LivePlots.PlotWidgetFactory(setup)
    Bases: object
```

The PlotWidgetFactory defines a simple interface for creating instances of previously defined LivePlotWidgets.

2.5 FAQ

Below known issues, their possible causes and subsequent fixes are listed. The fixes are ordered by decreasing likelihood so go from top to bottom retrying if the problem has been solved after every step.

2.5.1 Why don't I see real data?

Whenever the GUI is launched the setup tries to connect to all sensors. If that fails only simulated measurement values are shown. This can happen for the following reasons:

- The hardware does not have power.
 - 1. Check whether the experiment is plugged in.
 - 2. Check whether the power switch on the back is turned on.
 - 3. Check whether the fuses are intact.
- One of the USB devices is not connected.
 - 1. Connect the setup to your computer and validate that new devices are registered. Check the serials entry in Utility\config.json to see how many devices are expected to connect.
 - 2. When working on Windows, it can happen that two devices are registered under the same comport ID. To check, open the device manager and see if a comport ID appears twice. To fix either reassign one of the comport IDs manually or simply reboot your computer.

2.5.2 Why does the temperature difference decrease when heating?

- The temperature sensors are registered in the wrong order.
 - 1. Navigate to the dropdown menu on the top left Configuration -> Reverse Temperature Sensors.

2.5.3 Why does the flow controller not deliver sufficient flow?

- The pressurized air supply has only limited pressure.
 - 1. Open the pressure reduction valve a bit more until 100slm of flow can be delivered.

2.5.4 Why ... ?

- The program runs, the sensors are connected but nothing works as expected.
 - 1. Turn on debug mode, which allows you to view real time logs. See section *Debugging*.

PYTHON MODULE INDEX

u

Utility.Logger,8

INDEX

Symbols

```
_calibrate_temperature()
       (GUI.MainWindow.MainWindow
                                      method),
       15
_change_competition_mode()
       (GUI.MainWindow.MainWindow
                                      method),
       15
_convert_humidity() (Drivers.SHT.SHT method),
       13
_convert_temperature()
                              (Drivers.SHT.SHT
       method), 13
_go_to_next_view()
       (GUI.MainWindow.MainWindow
                                      method),
       15
_go_to_previous_view()
       (GUI.MainWindow.MainWindow
                                      method),
       15
_measure_normal_mode() (setup.Setup method), 6
_measure_simulation_mode()
                                   (setup.Setup
       method), 6
_reset_plots()
                   (GUI.MainWindow.MainWindow
       method), 15
_reset_temperature_calibration()
       (GUI.MainWindow.MainWindow
                                      method),
       15
_reverse_temperature_sensors()
       (GUI.MainWindow.MainWindow
                                      method),
       15
_save_measurement_buffer()
       (GUI.MainWindow.MainWindow
                                      method),
       15
_setup_measurement_buffer()
                                   (setup.Setup
       method), 6
_start_recording()
       (GUI.MainWindow.MainWindow
                                      method),
       15
_stop_recording()
       (GUI.MainWindow.MainWindow
                                      method),
       15
_toggle_massflow()
       (GUI.MainWindow.MainWindow
                                      method),
       15
```

_toggle_output()(<i>GUI.MainWindow.MainWindow</i>
method), 15
_toggle_setpoint()
(GUI.MainWindow.MainWindow method),
15
_update_lcds()(GUI.CustomWidgets.Widgets.StatusWidget
<i>method</i>), 16
_update_process_values()
(GUI. Custom Widgets. Widgets. Competition Widget
<i>method</i>), 16
_update_progress()
(GUI. Custom Widgets. Widgets. Competition Widget
method), 16

А

```
close() (setup.Setup method), 6
CompetitionWidget (class in
            GUI.CustomWidgets.Widgets), 16
connect() (Drivers.SFX5400.SFX5400 method), 14
connect() (Drivers.Shdlc_IO.ShdlcIoModule
            method), 11
connect() (Drivers.SHT.EKS method), 12
connect() (Drivers.SHT.SHT method), 13
connect_sensor() (Drivers.SHT.SHT method), 13
```

```
connect_sensors() (Drivers.SHT.EKS method), 12
```

D

DeviceIdentifier	(class	in
Drivers.DeviceIdentifier), 10		
disable_output()(<i>setup</i> .	Setup method), 6	

disconnect() (Drivers.SFX5400.SFX5400 method), Mode (class in setup). 8 14 disconnect() (Drivers.Shdlc IO.ShdlcIoModule method), 11 disconnect() (Drivers.SHT.EKS method), 12 disconnect() (Drivers.SHT.SHT method), 13

Е

EKS (class in Drivers.SHT), 12 enable_output() (setup.Setup method), 6

F

FancyPointCounter (class GUI.CustomWidgets.Widgets), 16

G

get analog input() (Drivers.Shdlc_IO.ShdlcIoModule method), 11 get analog output() (Drivers.Shdlc_IO.ShdlcIoModule method), 11 get_current_flow_value() (setup.Setup *method*), 6get_device_information() (Drivers.SFX5400.SFX5400 method), 14 get_digital_io() (Drivers.Shdlc_IO.ShdlcIoModule method), 11 (Drivers.Shdlc_IO.ShdlcIoModule get_pwm() method), 11

is_connected() (Drivers.SFX5400.SFX5400 method), 14 is_connected() (Drivers.Shdlc_IO.ShdlcIoModule method), 11 is_connected() (Drivers.SHT.EKS method), 12 is_connected() (Drivers.SHT.SHT method), 13

L

LivePlotSignal	(class				
GUI.CustomWidgets.LivePlots), 16					
LivePlotWidget	(class				
GUI.CustomWidgets.LivePlots), 17					
LivePlotWidgetCompe	etition (class				
GUI.CustomWidgets.LivePlots), 17					

Μ

MainWindow (class in GUI.MainWindow), 15 measure() (Drivers.SFX5400.SFX5400 method), 14 measure() (Drivers.SHT.EKS method), 12 measure() (Drivers.SHT.SHT method), 13 measure() (setup.Setup method), 6 MeasurementBuffer (class Utilin ity.MeasurementBuffer), 9

module Utility.Logger,8

Ο

open() (Drivers.DeviceIdentifier.DeviceIdentifier method), 10 open() (Drivers.PlatformBase.PlatformBase method), 10 open() (Drivers.SensorBase.SensorBase method), 12 open() (setup.Setup method), 6 Ρ

in

PlatformBase (class in Drivers.PlatformBase), 10 PlotWidgetFactory (class in GUI.CustomWidgets.LivePlots), 17

R

read_status_reg() (Drivers.SHT.SHT method), 13 RepeatTimer (class in Utility.Timer), 9 reset()(GUI.CustomWidgets.Widgets.CompetitionWidget method), 16 reset plot layout() (GUI.CustomWidgets.LivePlots.LivePlotWidget method), 17 reset_temperature_calibration() (setup.Setup method), 7 reverse_temp_sensors() (setup.Setup method), 7 run() (Utility.Timer.RepeatTimer method), 9

S

in

in

in

save measurement buffer() (setup.Setup method), 7 SensorBase (class in Drivers.SensorBase), 12 set_all_digital_io_off() (Drivers.Shdlc_IO.ShdlcIoModule method), 11 set analog output() (Drivers.Shdlc_IO.ShdlcIoModule method), 11 set_digital_io() (Drivers.Shdlc_IO.ShdlcIoModule method), 11 set_flow() (Drivers.SFX5400.SFX5400 method), 14 set_flow() (setup.Setup method), 7 set_kd() (setup.Setup method), 7 set_ki() (setup.Setup method), 7 set_kp() (setup.Setup method), 7 set_pid_parameters() (setup.Setup method), 7 set_pwm() (Drivers.Shdlc_IO.ShdlcIoModule method), 12 set_pwm() (setup.Setup method), 7 set_setpoint() (setup.Setup method), 7 set_temperature_calibration() (setup.Setup method), 7 Setup (class in setup), 6

```
module
                                           Util-
setup_custom_logger()
                            (in
        ity.Logger), 8
setup_status_bar()
        (GUI.MainWindow.MainWindow
                                        method),
        15
setup_tool_bar() (GUI.MainWindow.MainWindow
        method), 16
SFX5400 (class in Drivers.SFX5400), 14
ShdlcIoModule (class in Drivers.Shdlc_IO), 11
SHT (class in Drivers.SHT), 13
start_buffering() (setup.Setup method), 7
start_direct_power_setting()
                                    (setup.Setup
        method), 8
start_measurement_thread()
                                     (setup.Setup
        method), 8
start_pid_controller() (setup.Setup method), 8
StatusWidget
                            (class
                                             in
        GUI.CustomWidgets.Widgets), 16
stop_buffering() (setup.Setup method), 8
stop_measurement_thread()
                                     (setup.Setup
        method), 8
```

U

V